

# VISIONAIRE™ 5

5 LPM At-Home Oxygen Concentrator



## Maintenance-Free Design

With a unique design that eliminates preventative maintenance schedules, get the most of your stationary inventory without having to worry about extra patients visits and replacement filter costs.



### CAIRE INVESTS IN QUALITY

Design changes in November 2013 drastically increased overall durability, and improved the reliability of valves and sieve beds. This has reduced failure rates and has set a new standard for high-performing 5 LPM concentrators.



### REDUCED COST OF OWNERSHIP

Combining low failure rates with no preventative maintenance schedule creates a high-performing stationary with a low total cost of ownership.



### LIGHTWEIGHT AND POWER-EFFICIENT

At only 30 lbs (13.6 kg) with side and top handles, the VisionAire™ 5 is lightweight and easy to deliver. Its 290 W power draw at all settings is incredibly efficient.



#### CAIRE Contact Information:

Providers Toll Free  
1.800.482.2473  
Patients Toll Free  
1.877.704.0878  
[www.cairemedical.com](http://www.cairemedical.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.  
© 2020 CAIRE Inc. All Rights Reserved.

# VISIONAIRE™ 5

5 LPM At Home Oxygen Concentrator

## Includes

- Oxygen Concentrator
- Humidifier Adaptor Tube

## Product Features

- Maintenance-free design
- Lockable flow meter
- Lightweight and easy transport
- Dual-function hour meter
- Well-defined audible alarm system for easy diagnosis



## Specifications

**FLOW RATE**  
0.5 - 5.0 LPM

**DIMENSIONS**  
20.8 x 14.1 x 11.5 in  
(52.8 x 35.8 x 29.2 cm)

**WEIGHT**  
30 lbs (13.6 kg)

**SOUND LEVEL**  
45 dB(A) (Typical)

**O<sub>2</sub> CONCENTRATION**  
90 % (+5.5 / -3 %)

**OPERATING ENVIRONMENT**  
41 to 104°F (5 to 40°C),  
15-95% humidity

**STORAGE ENVIRONMENT**  
-13 to 158°F (-25 to 70°C)  
15-95% humidity

**POWER CONSUMPTION**  
290 W (All Flow Rates)

**ALARMS**  
Power Failure, High Pressure,  
Low Pressure, No Flow, Low O<sub>2</sub>  
Concentration (Optional)

**MAXIMUM TUBING**  
50 ft (15.2 m)

**MAINTENANCE SCHEDULE**  
None

**WARRANTY**  
3 Years

## Available Accessories



### HUMIDIFIER ADAPTOR & BOTTLE

Humidifier Adaptor (9 in/23 cm)  
Ask your doctor or therapist if  
humidification is right for you.



### NASAL CANNULAS

Various lengths are available.  
(7 ft and 25 ft)  
Talk to our Customer Service team  
to learn more.



### OXYGEN TUBING

Comes in: 7 ft (2.1 m), 25 ft (7.6 m),  
and 50 ft (15.2 m).



**CAIRE Contact Information:**  
Providers Toll Free  
1.800.482.2473  
Patients Toll Free  
1.877.704.0878  
[www.cairemedical.com](http://www.cairemedical.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.  
© 2020 CAIRE Inc. All Rights Reserved.

ML-CONC0082 D