



OXYGEN UNLIMITED

NUVO Series

Nuvo Standard

Nuvo 8

Nuvo 10



Nidek Medical's Nuvo Oxygen Concentrator carries forward our commitment of providing reliable, cost-effective and quality products.

An Oxygen Concentrator is a machine that generates medical oxygen. It runs on electricity and the only maintenance required is the regular cleaning of filters.

The Nuvo Oxygen Concentrator has been designed to provide a combination of enhanced features, reliability of a proven technology and unsurpassed quality.

The Nuvo Oxygen Concentrator is the largest selling concentrator in India. We promise you complete satisfaction and convenience.





NUVO Standard

- Quiet Operation - less than 48 dBA
- Lightweight - 50 lbs. - Easy to Transport and Maneuver
- Slender & Sleek Cabinet Design with Integrated Handles
- Quick Snap Rear Panel allows easy access to filter, hourmeter & battery
- Superior Grade of Molecular Sieve
- Patented Maintenance Free Rotary Poppet Valve
- Exclusive Moisture Reducing Wick Assembly
- Reinforced Base Substructure and Recessed Casters to guard against breakage
- Spring Mounted High Efficiency Thomas Nexus Compressor
- Patented RPSA Technology

Specifications

Size (LxWxH)	394x396x706 mm
Weight	22.7 Kgs
Average Power Consumption	420 Watts
Safety Alarms	Power Failure, Process Failure High & Low Pressure
Oxygen Purity	90% (+5%/-3%)
Output Pressure	7 psi
Oxygen Output	0 to 5 LPM
Filters	Cabinet, Compressor Intake & Bacteria
Biflow	Optional
Oxygen Monitoring System	Optional



NUVO 8

- Liter Flow - 1 to 8 LPM
- Ideal for High Flow Patients
- Low Maintenance Costs
- Made in the USA
- Built on the Reliable Nuvo Platform
- Easy Access to filter, battery, and hourmeter
- Superior Grade of Molecular Sieve
- Spring Mounted High Efficiency Thomas Nexus Compressor

Specifications

Size (LxWxH)	394x396x706 mm
Weight	25.2 Kgs
Average Power Consumption	500 Watts
Safety Alarms	Power Failure, Process Failure High & Low Pressure
Oxygen Purity	90% (+5% /-3%)
Output Pressure	15 psi
Oxygen Output	1 to 8 LPM
Filters	Cabinet, Compressor Intake & Bacteria
Biflow	Optional
Oxygen Monitoring System	Included



NUVO 10

- Liter Flow - 1 to 10 LPM
- Ideal for High Flow Patients
- Low Maintenance Costs
- Made in the USA
- Built on the Reliable Nuvo Platform
- Easy Access to filter, battery, and hourmeter
- Superior Grade of Molecular Sieve
- Spring Mounted High Efficiency Thomas Nexus Compressor

Specifications

Size (LxWxH)	394x396x706 mm
Weight	29.26 Kgs
Average Power Consumption	600 Watts
Safety Alarms	Power Failure, Process Failure, High & Low Pressure
Oxygen Purity	90% (+5% /-3%)
Output Pressure	15 psi
Oxygen Output	2 to 10 LPM
Filters	Cabinet, Compressor Intake & Bacteria
Biflow	Optional
Oxygen Monitoring System	Included



Nidek Medical India Private Limited

Head Office: 238B A.J.C. Bose Road, Kolkata 700 020
Tel: (033) 4017 7000 (30 Lines), Mob: 93395 66444

Delhi: 322, 3rd floor, Ansal Chamber II, 6, Bhikaji Cama Place, New Delhi 110 066
Tel: (011) 2616 5742, Mob: 93501 83222

Chennai: Plot No.1966, J- Block,,3rd Street, 13th Main Road, Anna Nagar (West), Chennai 600 040
Tel: 044 4266 3662

Mumbai: 211 B, Damji-Shamji Industrial Estate, 9 LBS Marg, Kurla (W), Mumbai 400 070
Tel: (022) 4215 5441 / 80809 14400

Hyderabad: H.NO-3-41-66, Sri Laxmi Nilayam, South Block, West Marredpally, Secunderabad 500 026
Tel: (040) 4015 0243, 93468 32266

Lucknow: 1/44 Vijayant Khand Gomti Nagar, Lucknow 226 010
Tel: (0522) 3191944

Chandigarh: 2nd Floor, SCO# 355, Sector 32-D, Chandigarh 160 031
Tel: 93562 78899

E-mail: sales@nidekindia.com, Website: www.nidekindia.com

