



10 litres

Product Code : MAOXY10
 (AE-10)

| SIZE (cm) | Length | Width | Height | kg NW | D M | PCS /CTN |
|-----------|--------|-------|--------|----------|------------|-------------|
| MAOXY10 | 37.2 | 34 | 61.2 | 24 kg | 44.5x40x68 | 1 |

Specifications :

| | |
|-----------------|--------------------------|
| Flow rate | 1-10L/MIN |
| Concentration | 90-96% |
| Power | 610W |
| Sound Level | <50db |
| Outlet Pressure | 100 Kpa +10% or 14.5 psi |
| Net Weight | 24 kg |
| Gross Weight | 26 kg |
| Warranty | 2 year |

- 2 Humidifer bottles, Dual flow
- Power failure Alarm; High Temperature Alarm and Voltage alarm
- Intelligent Self checking system: Alarm and Error Indication on Screen
- Long service life with stable performance.
- Suitable for 24 hours working
- Multi level filter & 2.5 micron Bacterial Filter
- Built-in Nebulizer
- Resetable Circuit breaker, Flame Breaker, Maintenance Alert.
- Most quiet machine. Best choice for indoor bedroom usage

2
 years
 Warranty*

CE 0123