

OXYGEN CONCENTRATOR

10 Litres



SINGLE FLOW

MCOXY10





Oxy-med[®] **eco**
O₂ for healthy living

OXYGEN CONCENTRATOR

10 Litres

SINGLE FLOW

SPECIFICATIONS

Flow rate	1-10L/MIN
Concentration	93+3%
Power	770W
Sound Level	≤55db
Outlet Pressure	70-80 Kpa
Operating Pressure	100-190 Kpa
Oxygen output Pressure	60kPa-80kPa
Power failure Alarm	Yes
Over-heat Alarm	Yes
Low Purity Alarm	Yes
Voltage Alarm	Yes
Oxygen Purity percentage display	Yes
Designed for 24 hour working	Yes
Net Weight	20 kg
Gross Weight	21.50 kg
Warranty	2 years