

# RESmart GII CPAP / Auto CPAP System

Resmart GII is a brand new intelligent CPAP device, suitable for home and clinic therapy of Obstructive Sleep Apnea (OSA). Its intelligent platform helps you breath comfortably and enjoy your sleep.





## RESmart GII CPAP / Auto CPAP System



#### 2.4 / 3.5-inch colorful LCD display

With an exquisite control knob, RESmart GII series bring you unsurpassed operational experience.

#### SpO<sub>2</sub> Kit

Real-time blood oxygen saturation and pulse rate display.

#### iCode QR

Compliance data is encoded by QR code. Search "BMCares" on Google Play or App Store.









	RESmart GII CPAP	RESmart GII Auto CPAP
Cover Color	Silver & Grey	Silver & Grey
LCD Screen	2.4 inches	2.4/3.5 inches
Mode	CPAP	CPAP/Auto
Pressure Range	4-20 cm H₂O	4-20 cm H₂O
Pressure Increment	0.5 cm H <sub>2</sub> O	0.5 cm H₂O
Reslex	Off/1/2/3	Off/1/2/3
Ramp	0-60 minutes	0-60 minutes
Sound	<30 dBA	<30 dBA
Alert	Power Fail/ High Leak	Power Fail/ High Leak
Humidity Level	Off/1-5	Off/1-5
Auto On/Off	Off (Default)	Off (Default)
Data Storage	SD Card	SD Card
Dimensions (L×W×H)	6.6" x 7.7" x 4.6" (170 × 196 × 118 mm) 11.4"×7.7"×5.3" (290 ×196 ×134 mm) with H60	
Weight	< 3.3 lb (< 1.5 kg) < 5.5 lb (< 2.5 kg) with H60	
Power Supply	Input range: 100 - 240 V AC, 50/ 60 Hz, 2.0 A max	

### BMC Medical Co., Ltd.