

OXYGEN CONCENTRATOR MODEL:JAY-5AW



OXYGEN CONCENTRATOR MODEL:JAY-5AW



USA imported humidifier,6PSI pressure release



Innovative "cooling" technology to improve reliability and life of the unit



Big LCD display
24H*365D work continuously



Specification

Flow Rate	0-5LPM
Purity	93%(±3%)
Outlet Pressure	0.04-0.07MPA(6-10PSI)
Sound Level	≤40db
Power Consumption	≤350w
LCD Display :	Switch times, operating pressure ,present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm, high&low pressure alarm
Product Weight	16 Kgs
Product Size	348x280x510mm



CORPORATE HEADQUATER

HEYER Medical AG

Carl-Heyer-Str. 1/3, 56130 Bad Ems, Germany
www.heyermedical.de

HEYER Medical AG

Carl-Heyer-Str. 1/3, 56130 Bad Ems, Germany
Tel: +49 2603 791 0
Fax: +49 2603 791 200
sales@heyarmed.com