

Bpure™ BiPAP Device

Advanced **Non-Invasive** Ventilation **Solutions**



Smarter Ventilation. Better Comfort.

BPure™ BiPAP Device is an advanced non-invasive ventilation system designed to provide effective respiratory support with intelligent pressure control, superior patient comfort, and reliable therapy performance.

Developed for hospitals, sleep labs, respiratory clinics, and homecare environments, BPure™ delivers precise bilevel ventilation therapy for patients with sleep-disordered breathing and respiratory insufficiency.

Advanced Features

- ✓ Auto-S, APAP, CPAP, BiPAP / S / T / ST/VAPS
- ✓ Intelligent Dual Pressure Ventilation
- ✓ Advanced Trigger Sensitivity Technology
- ✓ Real-Time Respiratory Monitoring
- ✓ Automatic Leak Compensation
- ✓ Heated Humidification Support
- ✓ Smart Backup Respiratory Rate
- ✓ Comfortable Pressure Relief System
- ✓ Auto Start & Stop Function
- ✓ Quiet & Stable Operation
- ✓ Modern Color LCD Interface

Intelligent Therapy Technology

BPure™ continuously monitors patient breathing patterns and automatically adjusts inspiratory and expiratory pressure levels to enhance therapy effectiveness and improve patient compliance. The advanced airflow system helps reduce therapy discomfort while maintaining stable respiratory support during sleep therapy.

Technical Specifications

Specification	Details
Product Name	Bpure™
Therapy Mode	Auto-S, APAP, CPAP, BiPAP / S / T / ST/VAPS
Maximum Pressure	30 cmH ₂ O
Humidity Levels	1-5 Levels
Pressure Relief	Adjustable
Display	4.3 inch Color LCD Screen
Leak Compensation	Yes
Auto Start / Stop	Yes
Humidification	Heated Humidifier
Application	Non-Invasive Ventilation
Size	215mm*167mm*145mm

Clinical Applications

- ✓ Obstructive Sleep Apnea (OSA)
- ✓ COPD Respiratory Support
- ✓ Hypoventilation Therapy
- ✓ Sleep-Disordered Breathing
- ✓ Homecare Ventilation
- ✓ Hospital Respiratory Management
- ✓ Long-Term Respiratory Therapy

Why Choose Bpure™?

Enhanced Patient Comfort

Advanced bi-level pressure delivery improves breathing comfort and therapy acceptance.

Intelligent Ventilation Support

Smart respiratory monitoring ensures stable and adaptive therapy.

Reliable Performance

Engineered for continuous respiratory support with consistent pressure delivery.

Ideal For

- ✓ Hospitals & Critical Care Units
- ✓ Sleep & Respiratory Clinics
- ✓ Homecare Providers
- ✓ Medical Equipment Distributors
- ✓ Respiratory Therapy Specialists

Better Breathing. Better Sleep. Better Life. BPure™ BiPAP Device



Ace Cathtech Pvt. Ltd.

Plot No 232, Sector 68, IMT Ballabhgarh, Faridabad Haryana -121004

Customer Care No.. +91 9266851166

Email: info@acecathtech.com

Website: www.acecathtech.com