



A comfortable mask improves adherence to treatment,
leading to faster well-being of the patients



AirFit™ N20



AirFit™ F20



Introducing **AirFit™ 20series**

Perfect Fit. Maximum Comfort.

AirFit™ N20



AirFit™ F20



- ✓ Adapts to your face for perfect fit
- ✓ Maximum seal even in high pressure & movements



ResMed Pty Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

Distributed by: ResMed India Pvt.Ltd.
205, 2nd Floor, Southern Park, Plot D2,
Okhla Canal, Saket, New Delhi 110017 India.
Tel: 011- 4781 3100, email: enquiry@resmed.co.in



A revolutionary technology

that allows the mask to adapt
to your face and provide maximum seal

Infinity Seal™ Cushion Technology
Designed to provide consistent support to ensure a confident seal regardless of facial type, therapy pressure or mask movement.

Soft nasal bridge
Shapes to a patient's nasal bridge to foster greater comfort and minimise the likelihood of red marks.

Plush silicone surface
Provides a softer cushion texture around the sensitive nasal region for a comfortable seal and easy fit.

Adaptive wings
Accommodate a wide range of facial types by curving to the profile of each wearer.

Supportive rim
Designed to deliver stability and security to the adaptive areas of the cushion for a robust seal.

Contoured chin region
Aims to add increased compliance to the chin area for greater comfort and seal performance during movement.

Polished silicone finish
Crafted to seal around the mouth and jawline for a smooth, reliable fit while preventing cushion ride up.



Supporting features for all-round comfort & superior experience

Visual freedom

Compact under eye frame, provides full visual freedom

Comfortable

Soft material that's durable and easy to use Velcro tabs

Soft frame

Combination of plastic frame + foam + fabric, provides built-in, no fuss comfort

High strength

Tested upto 400N (40kg) of force applied to frame arm attachments, delivers strength and flexibility.

Upgraded magnetic clips



Colour to colour usability guide with stronger grip

360° movement freedom

360° rotation at elbow-to-frame & elbow-to-hose connection

AirFit™ N20 & AirFit™ F20

fits **>95%** of all patients*



*Source: ResMed AirFit N20 Internal global fitting study of 159 existing ResMed patients, conducted 12/1/2015 : ID A3657023. ResMed AirFit F20 Internal International fitting study of 80 ResMed and non-ResMed patients, conducted between 08/03/2015 - 22/03/2015. Preliminary patient study - data on file: ID A3657701.