



Cpure™ CPAP Device

Smart Sleep Therapy Solutions



Advanced CPAP Therapy for Better Sleep

The CPAP Series is designed to deliver reliable, intelligent, and comfortable sleep therapy for patients suffering from Obstructive Sleep Apnea (OSA) and respiratory sleep disorders. Equipped with advanced pressure algorithms, humidification control, and smart monitoring functions, these devices provide high-performance respiratory care for hospitals, sleep labs and homecare providers.

Comfortable. Intelligent. Reliable.

CPure™ CPAP Device is designed to provide advanced sleep therapy with intelligent pressure management, enhanced patient comfort, and reliable respiratory support for Obstructive Sleep Apnea (OSA) therapy.

The sleek and modern CPure™ system combines smart monitoring technology with easy-to-use controls, making it ideal for hospitals, sleep clinics, distributors, and homecare patients.

Key Features

- ✓ CPAP / APAP / AUTO-S Therapy Modes
- ✓ Intelligent Pressure Adjustment
- ✓ Advanced Leak Compensation
- ✓ Heated Humidification System
- ✓ Automatic Start & Stop Function
- ✓ Comfortable Pressure Relief Technology
- ✓ Real-Time Respiratory Monitoring
- ✓ Modern Colour Display Interface
- ✓ Quiet & Stable Performance
- ✓ Compact and Elegant Design

Technical Specifications

Specification	Details
Product Name	Cpure™
Therapy Mode	CPAP / APAP / AUTO-S
Maximum Pressure	20 cmH ₂ O
Humidity Levels	1-5 Levels
Pressure Relief	1-3 Levels
Display	Color LCD Screen
Leak Compensation	Yes
Auto Start / Stop	Yes
Data Management	Yes
Application	Sleep Apnea Therapy

Clinical Applications

- ✓ Obstructive Sleep Apnea (OSA)
- ✓ Respiratory Sleep Disorders
- ✓ Homecare Respiratory Support
- ✓ Sleep Therapy Clinics
- ✓ Hospital Respiratory Care

Why Choose Cpure™?

Enhanced Patient Comfort

Soft airflow management and intelligent pressure adjustment help improve patient compliance and comfort during therapy.

Smart Monitoring Technology

Advanced respiratory event recognition and data tracking support efficient therapy management.

Reliable Respiratory Support

Designed for continuous therapy performance with stable pressure delivery and safety features.

Applications

- ✓ Obstructive Sleep Apnea (OSA)
- ✓ Respiratory Sleep Disorders
- ✓ Long-term respiratory therapy
- ✓ Homecare respiratory support
- ✓ Hospital and clinical respiratory management

Better Sleep. Better Breathing. Better Life.



ACE CATHTECH

Quality is key

Ace Cathtech Pvt. Ltd.

Plot No 232, Sector 68, IMT Ballabhgarh, Faridabad Haryana -121004

Customer Care No.. +91 9266851166

Email: info@acecathtech.com

Website: www.acecathtech.com