



The BestFit 2 full face mask

A vented full face mask designed for use with CPAP or BiPAP devices.

Adaptive
cushion

Flexible
frame

CPAP
compatible

BestFit 2 is part of ResMed India's BestFit & AcuCare mask range.

Designed around a simple, adaptable mask system

ResMed India lists the BestFit 2 with an adaptive cushion, flexible frame and CPAP compatibility.

Adaptive cushion
Part of the mask's sealing interface.



Flexible frame
Supports the mask assembly and is listed as a core feature by ResMed India.

Full face coverage
A full face mask interface.

Vented design
ResMed India describes BestFit 2 as a vented mask.

Three sizes
Small, Medium and Large are listed.

Main mask components

A visual overview of the components shown in the supplied BestFit 2 breakdown artwork.

ResMed BestFit 2 Full Face Mask

PART BREAKDOWN



Lightweight Design



Reliable Seal & Comfort



Quiet Exhalation



Easy to Fit & Adjust

BestFit 2 Full Face Mask
Designed to provide a secure seal and personalized fit for a comfortable therapy experience.





1 Headgear / Straps
Soft, adjustable fabric straps provide stability and hold the mask in place for a secure fit.

2 Adjustable Forehead Support
Provides additional stability and can be adjusted for a comfortable, customized fit.

3 Mask Frame
Sturdy frame that supports the cushion and connects all mask components.

4 Full Face Silicone Cushion
Soft, adaptive cushion that covers both the nose and mouth to deliver a comfortable and reliable seal.

5 Elbow / Swivel Connector
Connects the mask to the air tubing and allows freedom of movement.

6 Exhalation / Vent Ports
Specialized vent design allows exhaled air to escape quietly.

7 Headgear Clips
Secure the headgear to the mask frame. Easy to attach and detach.

To CPAP / BiPAP Tubing
Connects to the air tubing from your CPAP or BiPAP device.

Component	Function / role
Headgear / straps	Helps secure the mask on the head.
Forehead support	Provides upper-mask stability.
Mask frame	Supports the cushion and connects mask components.
Full face cushion	Provides the face-to-mask sealing interface.
Elbow / connector	Connects the mask to the therapy tubing.
Vent / exhaust area	Allows exhaled air to exit the mask system.

Sizing & product information

Choose the size that provides a stable and comfortable seal.

The BestFit 2 Full Face Mask is listed in Small, Medium and Large sizes. For fitting and sizing, use current ResMed product guidance or consult an authorised dealer.



Product codes

Size	Code	Dimensions
Small	80062	169 mm (H) × 106 mm (W) × 101 mm (D)
Medium	80061	169 mm (H) × 106 mm (W) × 101 mm (D)
Large	80060	169 mm (H) × 106 mm (W) × 101 mm (D)



BestFit 2

Full face mask
Vented mask for CPAP or BiPAP devices.

For fitting, replacement parts and therapy compatibility, refer to current ResMed documentation.

Source: ResMed India Best Fit 2 Full Face Mask product listing. Layout redesigned to follow the clean visual language of ResMed mask brochures.