



OXY365

5 Litre



Reliable Oxygen. Continuous Care.

The OXY 365 is a dependable, medical-grade oxygen concentrator engineered to deliver consistent and high-purity oxygen therapy for patients with respiratory needs. Designed with advanced technology, quiet performance, and user-friendly functionality, it is ideal for both homecare and clinical environments.





Integrated Atomization

Built in nebulization for direct medication delivery for respiratory care.



Continuous Operation

Stable oxygen therapy for longer use without interruption



Energy Efficient

Optimised power consumption for cost-effective operations.



High Oxygen Purity

Delivers medical grade oxygen with $93 \pm 3\%$ purity.



Advanced Benefits

Advanced Technology: Oil free compressor and using CECA molecular sieves for durability & efficiency. Designed to work for long hours without interruption.

Smart Safety Alarms: Alerts for power failure, low oxygen concentration and low pressure.

User-Friendly Interface: Simple control panels with Intelligent LED digital display.

Precise Atomization: Enables fine mist delivers with patient comfort.

Patented Technology: Backed by 8 EU-certified patents for reliability.

Low Noise Performance : Ensures patient comfort specially in nights.

One Purchase.
**Years of
Reliable Oxygen.**



Technical Specifications

Flow Rate Range	0-5 LPM adjustable
Oxygen Flow Rate	0.5-5 LPM
Oxygen Concentration	93-95 % $\pm$ 3%
Optimum Concentration	$\leq$ 1 minute
Oxygen Outlet Pressure	8 PSI $\pm$ 2
Noise Level	$\leq$ 45 dbA
Oxygen Display	Intelligent LED /Digital Display
Nebulization Port	Yes
Warranty	3 years or 15,000 Hours , Subject to Filter Replacement
Filters	Cabinet, HEPA Inlet , Bacteria
Altitude	Up to 10,000 Feet
Alarms	Power Failure, Low Concentration, Low Pressure
Operating Temperature	10-40 °C
Input Power	AC 220V, 50 Hz
Power Consumption	$\leq$ 350 W
Weight	$\leq$ 19 kg
Dimensions	340 *300*650 mm